


Diamond Range



Fibrelight
Laryngoscopes

Fibr Gr

Macintosh




Size 1 Infant 50530

Size 2 Child 50531

Size 3 Adult 50532

Size 4 Large 50533

Seward



Size 1 Infant 50538

Size 2 Child 50539

Miller



Size 0 Neonate 50534

Size 1 Infant/Child 50535

Size 2 Child/Adult 50536

Size 3 Large 50537

Robertshaw



Size 0 Neonate 50540

Size 1 Infant 50541

Fibrelight System

Lifelong Quality and Performance

Features

In 1972, Penlon designed the world's first Fibrelight Laryngoscope range. The **Diamond Laryngoscope Range** has been developed to meet the latest international fibrelight laryngoscope standard, ISO 7376-3, the Green Standard, to give a fully interchangeable, internationally compatible laryngoscope system. And the Diamond Range now includes the Fibrelight McCoy (see the McCoy brochure for full details).

- ◆ **Excellent fibre optics** - able to withstand repeated autoclave treatment
- ◆ **Five year guarantee** on light guide and lifetime guarantee on all blades and handles (excluding bulb)*
- ◆ **Removable light guide design** - extends blade life, and allows easier cleaning
- ◆ **Penlon** excellence in quality and design - sets the standards by which all others are judged
- ◆ **Quality** you can trust
- ◆ **Superior polished finish** - easy to clean and sterilise



Bulb Pack

Pack of Five 50542

Boxed Sets

| | |
|---|-------|
| Mil 3, Mil 2, Sew 1, Std Handle | 50374 |
| Mac 3, Mac 4, Sew 1, Std Handle | 50375 |
| Mac 1, 2, 3, 4 & Std Handle | 50376 |
| Special Custom sets available (Handles + up to 4 blades) | P.O.A |



Penlight 50526



Standard 50525



Stubby 50527

Technical Specification

| | | | | | |
|------------|-------|-------|--------|-------|-------|
| | | | | | |
| Mac 1 | x | y | Mil 0 | x | y |
| Mac 2 | 95 | 42 | Mil 1 | 76 | 35 |
| Mac 3 | 107 | 42 | Mil 2 | 102 | 35 |
| Mac 4 | 130 | 42 | Mil 3 | 152 | 37 |
| | 150 | 42 | | 192 | 37 |
| Macintosh | | | Miller | | |
| | | | | | |
| R/Shaw 0 | x | y | Sew 1 | x | y |
| R/Shaw 1 | 96 | 38 | Sew 2 | 110 | 37 |
| ----- | ----- | ----- | ----- | ----- | ----- |
| ----- | ----- | ----- | ----- | ----- | ----- |
| Robertshaw | | | Seward | | |

All dimensions are in mm and are approximate

Spares

| Light Guide Type / Size | 0 | 1 | 2 | 3 | 4 |
|-------------------------|-------|------|-------|-------|-------|
| Macintosh | ----- | 7265 | 7265 | 7267 | 7268 |
| Miller | 7269 | 7270 | 7271 | 7272 | ----- |
| Robertshaw | 7275 | 7276 | ----- | ----- | ----- |
| Seward | ----- | 7273 | 7274 | ----- | ----- |

| | |
|-----------------------------|-------|
| Click Screw | 7211 |
| Lamp Capsule Std / Penlight | 4410 |
| Lamp Capsule Stubby | 4411 |
| Bulb (Pack of Five) | 50542 |

**The Penlon Laryngoscope is guaranteed throughout its normal useful life against mechanical failure due to defects in design, materials or workmanship. The guarantee does not apply to lamps. Fibrelight guides are guaranteed against total light failure within a period of five years. However note that after many sterilisations the light output of the fibre bundle may decrease. Contact the Customer Service Department at Penlon Ltd for details of our repolishing service.*



Penlon Ltd, Radley Road, Abingdon, Oxon OX14 3PH England
 Tel: +44 (0) 1235 547001, Fax: +44 (0)1235 547021,
 E-Mail: export@penlon.co.uk, uksales@penlon.co.uk
 Website: <http://www.penlon.com>