


# Maloney & Hurst Bougie

## For Esophageal Dilation



Rüsch bougies, for esophageal dilation, are available in both Hurst and Maloney styles. Designed with silicone casings for easy insertion when wetted with sterile water. The dry areas of the casing remain tacky to allow for a sure grip. Weighted with Tungsten for easy disposal when necessary. The end cap has a permanent french size marking for easy identification.

### Specifications:

- ◇ Tungsten filled (Mercury-Free)
- ◇ Depth Markings
- ◇ Radiopaque throughout entire length
- ◇ Easy Disposal
- ◇ Clear marking along the side of the bougie

### Benefits:

- ◇ Silicone material becomes slippery when wetted with sterile water, no addition lubrication required.
- ◇ Maloney taper is gradual and smooth
- ◇ Convenient carrying case available
- ◇ Available in a full range of sizes

### BOUGIES



| Maloney | Size | Hurst  |
|---------|------|--------|
| 215020  | 20   | 215520 |
| 215022  | 22   | 215522 |
| 215024  | 24   | 215524 |
| 215026  | 26   | 215526 |
| 215028  | 28   | 215528 |
| 215030  | 30   | 215530 |
| 215032  | 32   | 215532 |
| 215034  | 34   | 215534 |
| 215038  | 38   | 215538 |
| 215040  | 40   | 215540 |
| 215042  | 42   | 215542 |
| 215044  | 44   | 215544 |
| 215046  | 46   | 215546 |
| 215048  | 48   | 215548 |
| 215050  | 50   | 215550 |
| 215052  | 52   | 215552 |
| 215054  | 54   | 215554 |
| 215056  | 56   | 215556 |
| 215058  | 58   | 215558 |
| 215060  | 60   | 215560 |

**Maloney Set** 215000

**Hurst Set** 215500

**Maloney Case** 215005

**Hurst Case** 215505



**Teleflex**  
MEDICAL

To place an order contact  
Customer Service at  
**1-800-553-5214**

ANESTHESIA  
2917 Weck Drive  
Research Triangle Park, NC ♦ 27709  
www.rusch.com  
P: 800 553 5214 ♦ F: 770 623 1829