

SALTER LABS®



Improved

Humidifier

#7100 and #7600 350 cc
Bubble Humidifier - Dry

ALL NEW 360° Rectangular Diffuser ports provide quiet operation with low system back pressure. New dome lid design incorporates refinements, resulting in high performance with rugged durability. All Salter bubble humidifiers provide, dependable, cost-effective performance. Ideal for use with all oxygen sources: concentrators, liquid, tanks, wall outlets...**Available with 3 PSI or 6 PSI Pop-off Alarm.**



Features

- NEW 360° molded, low resistance, high output diffuser head design
- NEW connecting tower design
- NEW large, quad wing nut connector
- NEW, audible pop-off alarm
3PSI or 6PSI
- Jar permanently etched, large print minimum/maximum water level markings
- Easy seal, translucent, double helix thread jar

Benefits

- Quiet, efficient hydration with less rain out. Resists clogging due to mineral build up from water.
- Greater resistance to damage from incidental impact
- Easy to grip, reduces chance of "cross-threading"
- Reduces the likelihood of concentrator damage. Audible warning of downstream occlusion.
- Easy to read, durable will not rub, wear or wash off
- Simple, sure, half turn on/off with positive seal. No gasket is required

Variety of uses... Homecare, Transport, Hospital, Extended Care Facilities, Hospice



The Salter Humidifier #7100 3psi pop off #7600 6psi pop off

Disposable Humidifiers from Salter Labs are "truly single patient use products designed and engineered for long term durability and complete patient satisfaction".



Start-up Kit
Available with 7100
or 7600 humidifier

New, large, easy to grip quad style wing nut. Resists cross threading, simplifies secure connection to oxygen supply source.


New, audible safety pressure pop-off alarm. Functions at low flow rates (3PSI or 6PSI).

New, easy-flow dome lid design combined with unique 360° diffuser helps to reduce rain-out.

Easy jar to lid assembly with half turn positive seal.

Stem and diffuser resists clogging due to mineral build-up on and inside diffuser

New, permanent, maximum-minimum water level indicators. Will not wash, wear or rub off.



Unique 360° rectangular diffusion ports operate quietly with improved hydration and less system back pressure.

Ordering Information

Humidifier	Units (per case)		Part No.	
	3PSI	6PSI		
350 cc bubble humidifiers 3PSI (brown lid)	50		7100	
350 cc bubble humidifiers 6PSI (black lid)	50			7600

For high flows (6-15 lpm) refer to our model #7900 Bubble Humidifier.

Start-up Kits

2 #1606 nasal cannulas, 1 #7100 or #7600 humidifier and 1 #2002 supply tube:

Kit with a 7' tube	20	7107	7607
Kit with a 25' supply tube	20	7125	7625
Kit with a 50' supply tube	20	7150	7650

Start-up Kits with smooth bore tubing available on request



100 W. Sycamore Road - Arvin, CA 93203
Phone: 661-854-3166
1-800-421-0024 (CA only)
1-800-235-4203 (U.S. & Canada)
Fax: 661-854-3850
Toll Free Fax: 1-800-628-4690 (U.S. & Canada)
International Fax: 661-854-3850
E-mail: salterlabs@us.salterlabs.com
www.salterlabs.com

CE
0482

ISO 9001/EN 46001 & ISO 13485