

NEW FROM WELCH ALLYN...

TWICE THE POWER, HALF THE WEIGHT



INTRODUCING THE WELCH ALLYN POWER HANDLE ...SMARTER TECHNOLOGY

For the first time in a power handle, Welch Allyn introduces smarter Lithium Ion technology to provide twice the power with just half the weight. The 719 Series Power Handle's new Lithium Ion technology powers any 3.5v instrument head for 110 minutes – twice as long as the conventional handle. And at half the weight and with an updated, ergonomic design, the 719 Series is more comfortable for any size hand. An indicator light warns you when the battery is running low, and the built-in AC adapter makes recharging simple and convenient. Made of durable CYCOLOY® construction, the 719 Series Power Handle will provide years of reliable quality.

WelchAllyn®
719 Series
Power Handle

LITHIUM ION TECHNOLOGY PROVIDES TWICE THE POWER WITH JUST HALF THE WEIGHT.

719 Series Power Handle Features

- Half the weight of current 3.5v handles
- Average 110 minutes of on-time between charging
- Smallest diameter of any 3.5v power handle
- Made from CYCOLOY® resin from GE plastics
- Low-battery indicator
- Lithium Ion battery technology
- ETL listed to UL 2601 and CAN/CSA C22.2 No. 601.1 medical electrical equipment standards
- Module section can be removed during use

Lithium Ion Battery Technology

- No “memory” effect
- Low battery indicator
- No need to discharge between charging cycles



(Handle shown at actual size)

AC Charging Module

- Easily removable from handle
- Retractable AC prongs
- Can be attached or detached to handle during exams

*CYCOLOY is a registered trademark of General Electric Company



Charger module can be detached to provide the smallest 3.5v handle available.



Charger module can be left attached for enhanced convenience and comfort for larger hands.



Retractable module can be plugged into any 110v outlet.

ORDERING INFORMATION

- 71900** 719 Series Power Handle, includes complete removable AC adapter
- 71950** 719 AC Charging Module
- 71960** 719 Replacement battery

U.S. and International patents pending

For more information about the Welch Allyn 719 Series Power Handle, contact Welch Allyn at **888-359-2394**.

WelchAllyn®
719 Series
Power Handle

4341 State Street Road, Skaneateles Falls, NY 13153-0220 Tel: (888) 359-2394, (315) 685-2679 Fax: (315) 685-2174
www.welchallyn.com infomedical@mail.welchallyn.com