

ZOLL

Total Protection for the First and Only Full-Rescue AED

Unbreakable, watertight, dustproof, chemical resistant and corrosion proof

Providing an excellent shock-resistant carrier for your ZOLL AED Plus™, the Pelican Case and its molded high-density foam provide total protection from the rigors of impact, vibration, and shock. Made of high-impact structural copolymer, it is extremely strong and durable. With its exclusive 0.25" (6.4 mm) neoprene O-ring, and ABS latches, this case seals perfectly. The Pelican Case also includes an automatic purge valve for quick equalization after changes in atmospheric pressure. The AED Plus Pelican Case is NATO-codified and tested to MIL C-4150J (Military Standard) and IP67 (water and dust ingress protection).

You break it, we replace it... forever.



Small Case
(19" x 15.4" x 7.6")
(48.3 cm x 39.1 cm x 19.3 cm)

AED+PLUS[®]

PELICAN CASE

ZOLL
Advancing Resuscitation. Today.™

**Waterproof
Dustproof
Crushproof**

**Total Protection for
the First and Only
Full-Rescue AED**

ZOLL offers two cases to meet your needs: a large case for easy storage of *CPR-D•padz* and *pedi•padz II* or a smaller case for the AED Plus alone. Each case comes complete with custom-molded foam inserts for a perfect fit.

Large (8000-0837-01)

Large Case with cutouts for AED Plus, *CPR-D•padz™*, or *pedi•padz® II*


Small (8000-0836-01)

Small Case with cutouts for AED Plus only

Inside Dimensions	17.9" x 12.8" x 6.8" (45.5 cm x 31.8 cm x 17.3 cm)	14.7" x 10.4" x 6" (37.3 cm x 26.4 cm x 15.2 cm)
Outside Dimensions	19" x 15.4" x 7.6" (48.3 cm x 39.1 cm x 19.3 cm)	16" x 13" x 6.9" (40.6 cm x 33 cm x 17.5 cm)
Temperature Ratings	Min: -10° F (-23° C) Max: 210° F (99° C)	Min: -10° F (-23° C) Max: 210° F (99° C)
Buoyancy	Floats in salt water with 40 lbs. (18.1 kg) load	Floats in salt water with 20 lbs. (9.07 kg) load
Color	Black	Black
Approvals	IP67 Military Standard C-4150J	IP67 Military Standard C-4150J



**A Rugged Case for
the First and Only
Full-Rescue AED**



ZOLL's AED Plus Large Pelican Case allows you to store your AED Plus and an extra set of *CPR-D•padz* with the option of also storing a set of *pedi•padz II*

ZOLL Medical Corporation
Worldwide Headquarters
269 Mill Road
Chelmsford, MA 01824
978-421-9655

For worldwide locations, visit:
www.zoll.com/contacts

Specifications subject to change without notice.
© 2006 ZOLL Medical Corporation. All rights reserved.
"Advancing Resuscitation. Today.", AED Plus and *CPR-D•padz* are trademarks of ZOLL Medical Corporation.
pedi•padz and ZOLL are registered trademarks of ZOLL Medical Corporation.
ZOLL Medical Corporation 269 Mill Road Chelmsford, MA 01824-4105 USA
9656-0173 Printed in USA 5M 7/06

